

ENDOLASER 120



HIGH PERFORMANCE THERAPEUTIC LASER

The 905 nm GaAs laser probe delivers pulses in nanoseconds, producing average powers of 100 mW and peak powers up to 100 W per diode. This results in a higher concentration of light energy (I_0), or photon density at tissue depth, without the risk of burning tissue.

PROVEN EFFECTIVENESS

The clinical application of light – produced by laser diodes with a power output up to 1600 mW – is scientifically well documented. The 905 nm (GaAs) pulsed laser is the most versatile for deep lying structures (joints, vertebrae, back (muscles) whereas the 808 nm (GaAlAs) is more appropriate for wound healing and tissue repair. Scientific evidence continues to demonstrate that pulsed light does have biological and clinical effects that are different from those of continuous wave (CW) light. Several studies revealed that the LLLT in pulsed wave mode of operation can better penetrate through the melanin and other skin barriers, supporting the hypotheses that pulsing is beneficial in reaching deep target tissue and organs [1].

EVIDENCE BASED CLINICAL GUIDELINES

The Endolaser provides the user with dosage recommendations for the effective treatment of pain and various musculoskeletal disorders.

MULTIPLE LASER PROBES AVAILABLE

For the treatment of smaller surfaces probes, there is a choice of single laser probe with a total power up of 100mW and 500mW are available. For the treatment of larger areas the cluster probe, combining 4 diodes with a total up to 1600 mW, can be used. The Endolaser automatically recognizes the type of probe that is connected.

QUICK ACCESS POWER BUTTONS

The output power of the Endolaser 120 can be easily set at different levels for easy dosage settings.

TARGET LIGHT

When treating a patient without probe-skin contact, the red target light helps the therapist to focus on the area to treat. The target light contributes to a safe treatment by indicating the irradiation area of the laser beam.

ERGONOMIC PROBES

The ergonomic probes can be used with or without skin contact. Each probe is equipped with a start-pause button for easy control. The indication light on top of the probe indicate the emission state (laser running or laser-ready).

ENDOLASER 120

1633901 ENDOLASER 120

with single probe holder



STANDARD ACCESSORIES

- Endolaser 120 base unit
- 2x Laser Safety Eyewear (goggles)
- Infosheet
- User manual (on CD)
- 3444290 Mains cable 230V-EUR
- 1x Probe holder (premounted)
- 1x Doorswitch cable

1633902 ENDOLASER 120

with double probe holder



STANDARD ACCESSORIES

- Endolaser 120 base unit
- 2x Laser Safety Eyewear (goggles)
- Infosheet
- User manual (on CD)
- 3444290 Mains cable 230V-EUR
- 2x Probe holder (premounted)
- 2x Doorswitch cable
- Interlink cable (premounted)

OPTIONAL ACCESSORIES

- | | | |
|---------|---------------------------------------|---------------------------------------------------------|
| 1632801 | LP100P | Laser probe 100 mW, super pulsed laser diode (905 nm) |
| 1632802 | LP500C | Laser probe 500 mW, continuous laser diode (808 nm) |
| 1632803 | CP4X100C | Cluster probe 4x100 mW, continuous laser diode (808 nm) |
| 1632804 | CP4X400C | Cluster probe 4x400 mW, continuous laser diode (808 nm) |
| 3444820 | Laser Safety Eyewear (set of goggles) | |

OUTPUT

- 100 mW
- 500 mW
- 4 x 100 mW
- 4 x 400 mW

PEAKPOWER

- 100 W



1632801



1632802



1632803



1632804



3444820

TECHNICAL SPECIFICATIONS

General

- | | |
|----------------------------------------------|-----------------|
| Mains voltage: | 100 - 240 Volt |
| Frequency: | 50/60 Hz |
| Max. power output: | 20 VA |
| Dimensions device (length x width x height): | 22 x 16 x 14 cm |
| Weight unit: | 1.6 kilograms |

- | | |
|-----------------------------------|---------------------------------------|
| Type of laser diode: | GaAs (pulsed) and GaAlAs (continuous) |
| Preprogrammed clinical protocols: | 22 |
| Number of favorites to save: | 20 |

Safety and performance standards
Medical device classification
93/42/

with

IEC 60601-1

Safety class according to IEC 60601-1
Applied parts
Classification of laser according to IEC60825-1 3B

IIa; Rule 9 Annex IX of EEC
This equipment complies all requirements of the Medical Device Directive (93/42/EEC).
General requirements for the safety of electrical medical systems.
Electrical safety class II
Type B Applied part.



ENRAF-NONIUS B.V.
Vareseweg 127
3047 AT Rotterdam
Nederland
T: +31-(0)10 - 203 06 00
E: info-nl@enraf-nonius.nl
www.enraf-nonius.nl

ENDOLASER 120

Powerful solution to eliminate pain,
reduce inflammation and
accelerate tissue healing

Fast & easy operation:
"with the **speed** of light"




Safe



Evidence Based Clinical **Guidelines**
incorporated


Drug-free and non-invasive

Effective




Smart laser monitoring system
-Long Life Time Laser Technology-
to ensure a **stable output**.
Now and in de the future

ENDOLASER 120



Automatic recognition of connected
laser probes



High Quality
beam characteristics




Full colour **touchscreen**

Automatic dosage / time adjustment
based on output power level



Ergonomically designed probes



Variety of probes:
CW and Pulsed



Easy to **clean** probes

Compact & light weight