



RT 55 TRANSVERSAL TENDER Physiotherapy Catalogue

Mediotronics Physical Medicine (Pty) Ltd

Awarded Tender items List

1. Heat Gun - RT55-A-16
2. Goniometer Set (Plastic) - RT55-A-18
3. Carry case for Goniometer Set - RT55-A1-18
4. Finger Goniometer (Stainless Steel) - RT55-A-19
5. Goniometer 180° 15-20cm (plastic) - RT55-A-20
6. Goniometer 180° 30cm (plastic) - RT55-A-22
7. Mobile Posture Mirror (Steel) - RT55-A-23
8. Electric Standing Frame, Adjustable (patient lifting) - RT55-A-29
9. Upper Extremity Workstation - RT55-A-34
10. Hydraulic Hand Dynamometer - RT55-A-37
11. 5 Section Height Adjustable Plinth (electronic) - RT55-A-54
12. 5 Section Height Adjustable Plinth (electronic) - RT55-A-55
13. Neurological Bobath Plinth (electronic) - RT55-A-57
14. Corner Steps / Stairs, with handrails - RT55-A-64
15. Active / Passive Motion Trainer, upper and lower body- RT55-A-65



Heat Gun **R1770** incl VAT & delivery RT55-A-16

2000W output

Carry Case with 4 nozzles

Reworking Spatula




6pc Goniometer Set **R650** incl VAT & delivery *excludes soft ca case RT55-A-18

Includes 15cm 20cm and 30cm 360° goniometers, finger goniometer,


SOFT CARRY CASE – R365 incl VAT - RT55-A1-18



Stainless Steel Finger Goniometer **R1200** incl VAT & delivery RT55-A-19



15cm Goniometer 180° **R180** incl VAT & delivery RT55-A-20



30cm Goniometer 180° **R245** incl VAT & delivery RT55-A-22



Mobile Posture Mirror **R3550** incl VAT & delivery RT55-A-23

0-360° Vertical Tilt, powder coated



Electric Standing Frame **R97850** incl VAT & delivery RT55-A-29

AkcesMed Lifter™ Electric Standing frame with lifting capabilities. Allows the patient to independently use the device by means of the remote control

Maximum weight capacity 120Kg



Upper Extremity Workstation **R32950** incl VAT & delivery RT55-A-34

Complete with ALL 11 activities, including, Bilateral Sander, extend-your-reach, E-Z Exercise board, Finger Extension, Peg Board, Shoulder Abduction Ladder, Shoulder Exercise Ladder, Manipulation board, Vertical ring Tree, Horizontal ring Tree and Stacking Cones

Easily Adjustable side trays, lockable cabinet and braked castors

Hydraulic Hand Dynamometer **R5200** incl VAT & delivery RT55-A-37

Saehan Hydraulic Hand dynamometer, in Carry case

0-90Kg force measurement, Adjustable hand with 5 positions for any hand size





5-Section Exam Couch **R25150** incl VAT & delivery

RT55-A-54

Electrically Adjustable height, and Central Lumbar Flexion via the foot pedal, Additional adjustments for Head Section, Leg Section and Arm supports

Maximum weight capacity 250Kg

Working Height from 46-91cm



3-Section Exam Couch **R23950** incl VAT & delivery

RT55-A-55

Electrically Adjustable height, and Central Lumbar Flexion via the foot pedal, Additional adjustments for Head Section and Leg Section, 198 x 70 x 91cm

Maximum weight capacity 250Kg

Working Height from 46-91cm



2-Section Bobath Plinth **R25950** incl VAT & delivery

RT55-A-57

Electrically Adjustable height by foot control with Additional adjustments for Head Section, 208 x 120 x 91cm

Maximum weight capacity 330Kg

Working Height from 46-91cm



Corner Steps **R15983** incl VAT & delivery

RT55-A-64

High Quality Lacquered Wood, With Sturdy Hand rails

Carpet Fixed on each providing better grip



Active / Passive Motion Trainer **R132980** incl VAT

delivery

RT55-A-65

Reck Motomed Loop Professional

Biofeedback and motivation

The MOTOMed Loop provides a detailed (biofeedback) during the training and it enables you to train even more effectively. The exciting motivational training games increase the training success and are fun to play.

MOTO med loop.la pro includes the following special features:

- Plastic-coated safety foot shells for easy disinfection / cleaning
- Leg guides with plastic-coated calf shells and disinfectable fastening straps
- Quick Fix « foot fastening system
- Plastic-coated handles with quick release system
- Gas spring support for height adjustment
- Pedal radius quick adjustment

Therapy options

- passive training (with motor)
- assistive training (motor-assisted)
- active training (with your own strength)
- Symmetry Training
- Coordination training (4-Segment)

Software

- big, user-friendly operating panel; color code system and big, easy-feel buttons
- language selection between 27 different languages
- MOTO max and TRAMPOLINE max motivation program
- therapy programs
- detailed feedback in real time
- post training analysis