


Spirobank II Plus



Available with Oximetry option

Spirobank II Plus is the latest generation of **portable spirometers** with **oximetry** option from MIR's Plus Line.

The integrated **Tag Reader** automatically identifies the **Flowmir Plus flow sensor**, new and Made in Italy turbine enhancing **low flow measurement** to ensure the highest standards of **accuracy** and **unmatched reliability**.

Its **320x240 color touchscreen** allows for faster and easier data entry and a more immersive user experience.

With its **high-definition screen**, Spirobank II Plus ensures clear and detailed views during every stage of the test, enhancing the interaction with the device.

Via **Bluetooth** or **USB**, Spirobank II Plus allows for the import of **Work Lists** from **MIR Spiro Software** and provides test results, significantly improving the overall efficiency of the diagnostic process.

Spirobank II Plus is designed to seamlessly integrate into modern clinical practices, offering **versatility** and **reliability** in every situation, from **Stand-Alone** use to connection to **Tablets** and **PCs**.

CHOOSE TEST

Massimiliano Arcangeli
AGE 33 | 190cm | 90kg | Male

Spiro **Oxi** **Data** (0 Test)

PRE PHASE

SPIROMETRY RESULT

Param.	Pred.	PRE#1	%Pred.	PRE#2	PRE#3
*FVC	4.77	1.61	34	1.61	1.61
*FEV1	4.10	1.54	38	1.54	1.54
*PEF	8.85	1.73	20	1.73	1.73
FVC	4.77	1.61	34	1.34	1.25

OXIMETRY RESULT

%SpO2	BPM
98,6	86,4
MIN 98 MAX 99	MIN 98 MAX 99

Best Features

- **Tag Reader integrated**
Reading and identification of data from the Plus Turbine RFID




- **Touchscreen Display**
320 x 240 Color LCD

- **FlowMIR Plus**
 - Automatic turbine type recognition
 - Factory calibration stored on RFID
 - Turbine traceability

- **Hardware**
New STM32 technology





- **Work List**
Bidirectional communication between the device and EMR/ MIR Spiro Software




● Spirometry

● Oximetry



Param.	Pred.	PRE#1	%Pred.	PRE#2	PRE#3
*FVC	4.77	1.61	34	1.61	1.61
*FEV1	4.10	1.54	38	1.54	1.54
*PEF	8.85	1.73	20	1.73	1.73
FVC	4.77	1.61	34	1.34	1.25



Mode of Use

- Spirobank II Plus has indeed **3 Modes of Use** in a single Device.

Spirobank II Plus is designed to integrate seamlessly into modern clinical practices.

It offers **versatility** and **reliability** in every situation, from Stand-Alone to use with tablets and PCs.



Stand Alone



With Tablet Connection



With PC Connection